

**Quantifying body
composition from CT
scans using commercially
available software:**

Applications in radiology




Authors

Varun Danda¹, Jessa Tunacao², Omid Shafaat², Clifford Weiss²

¹Sidney Kimmel Medical College, Philadelphia, PA, ²Johns Hopkins Hospital, Baltimore, MD

No relevant conflicts of interest to disclose.





Introduction

- Interest in body composition measurement techniques have skyrocketed over the last few years
- Measurement of this data has value in answering important research questions in many patient populations
 - Cancer patients
 - Transplant candidates and recipients
 - Other procedures (TIPS, EVAR)




Purpose

- Here we will provide an overview of a precise method to obtain body composition data from CT scans using the commercially available software OsiriX (Pixemo, Switzerland)




Methods

- OsiriX can be linked to an institution's DICOM network for easy retrieval of desired studies
- Patient data can be anonymized through the software to remain compliant with research protocols





Methods

- Once a study is opened, localize the L3 vertebral level using the sagittal view and locate a slice at the L3 vertebral level
- Proceed with segmentation using the axial view




Methods

- On the selected axial slice, use the closed polygon tool to draw a boundary between the abdominal muscle layer and the visceral adipose tissue layer as shown




Methods

- Next the “Grow Region 2D/3D Segmentation” tool is used to create ROIs with chosen upper and lower HU values
- The three tissue types chosen for segmentation with their HU thresholds reported in the literature are below:
 - Subcutaneous Adipose Tissue (SAT): **-190 to -30 HU**
 - Visceral Adipose Tissue (VAT): **-150 to -50 HU**
 - Skeletal Muscle (MM): **-29 to +150 HU**



Methods

- When segmenting skeletal muscle and SAT, select the “Out” option to ensure that only pixels outside the polygon are chosen to represent this tissue
- Similarly, segmentation of VAT requires selection of the “In” option
- Label each ROI appropriately




Results

- Preview Result when clicking
- Limit Results to selected ROI:
 - In Out
- Directly generate ROI when clicking
- Merge with existing Brush ROIs
- Propagate result in 4D
- Create ROI(s) in the original series
 - Brush ROI Polygon
 - # of points: Low High
 - ROI name:
- Generate a new series with:
 - Set Inside Pixels to:
 - Set Outside Pixels to:

Compute

Methods



Green: Subcutaneous Adipose Tissue (SAT)


Red: Skeletal Muscle (MM)

Yellow: Visceral Adipose Tissue (VAT)

- Example of final segmentation
- Data can be exported using the “ExportROIs” plugin which creates a .csv file
- Relevant data includes the total area of the ROIs by label as well as the HU distribution within the tissues


Reliability Test

- An attending radiologist, a PGY-1, and a MS-2 segmented and quantified VAT, SAT, and skeletal muscle at the L3 vertebral level from 140 CT scans using this method.
- The area data of all three tissues was exported from OsiriX for statistical analysis.
- An intra-class correlation coefficient (ICC) calculation demonstrated an ICC of **0.998** for SAT area, **0.999** for VAT area, and **0.995** for skeletal muscle area.




HU Density Plots

- HU data from skeletal muscle can be used to create density plots using the software R (R Core Team, Austria)
- Changes in muscle density can be measured over time and graphically visualized using these plots



Conclusions

- The method described appears to accurately measure body composition metrics from CT scan regardless of the training level of the reader.
- This technique can be useful in prognostic studies where body composition can be correlated to patient outcomes i.e. oncology or surgical patients.
- Studies that follow patients with serial CT scans can use this technique to detect longitudinal changes in body composition.



References

1. Aubrey, J., et al. "Measurement of Skeletal Muscle Radiation Attenuation and Basis of Its Biological Variation." *Acta Physiologica*, vol. 210, no. 3, 2014, pp. 489–497., doi:10.1111/apha.12224.
2. Mourtzakis, Marina, et al. "A Practical and Precise Approach to Quantification of Body Composition in Cancer Patients Using Computed Tomography Images Acquired during Routine Care." *Applied Physiology, Nutrition, and Metabolism*, vol. 33, no. 5, 2008, pp. 997–1006., doi:10.1139/h08-075.
3. Nattenmüller, Johanna, et al. "Prognostic Impact of CT-Quantified Muscle and Fat Distribution before and after First-Line-Chemotherapy in Lung Cancer Patients." *Plos One*, vol. 12, no. 1, 2017, doi:10.1371/journal.pone.0169136.
4. Pezaro, C, et al. "Sarcopenia and Change in Body Composition Following Maximal Androgen Suppression with Abiraterone in Men with Castration-Resistant Prostate Cancer." *British Journal of Cancer*, vol. 109, no. 2, 2013, pp. 325–331., doi:10.1038/bjc.2013.340.

